



#### CE CONICAL CENTRIFUGE TEST TUBES WITH SILKSCREEN PRINTED GRADUATION

**Material: Polypropylene**

These new lines with capacity of 15 ml and 50 ml, are CE marked conforming standard 98/79/CE, they are non-pyrogenic and autoclavable. Conceived for centrifuge, they are also suitable for biological sampling collection and transport.

#### Blue Silkscreen Printed Graduations:

Tube Capacity	Grad. Range	Sub-divisions
15 ml	2 - 15 ml	0.5 ml
50 ml	5 - 50 ml	2.5 ml

- Leakproof blue screw cap, knurling walls to make opening and closing easier. Conical bottom to allow complete recovery of sample liquid.
- Temperature resistant -80°C to +121°C

Art.	capacity	gamma radiation sterilised	packaging	quantity	Description		
CE CONICAL CENTRIFUGE TEST TUBES WITH SILKSCREEN PRINTED GRADUATION	<b>84006</b>	15 ml	n/a	Bulk	50	Blue Grad. Conical Test Tubes, 15 ml	<b>NEW</b>
	<b>84007</b>	15 ml	Sterile	Bulk	50	Blue Grad. Conical Test Tubes, 15 ml	<b>NEW</b>
	<b>84008</b>	15 ml	Sterile	Indiv. wrap	1	Blue Grad. Conical Test Tubes, 15 ml	<b>NEW</b>
	<b>84009</b>	15 ml	n/a	PS foam rack	50	Blue Grad. Conical Test Tubes, 15 ml	<b>NEW</b>
	<b>84010</b>	15 ml	Sterile	PS foam rack	50	Blue Grad. Conical Test Tubes, 15 ml	<b>NEW</b>
	<b>84011</b>	50 ml	n/a	Bulk	25	Blue Grad. Conical Test Tubes, 50 ml	<b>NEW</b>
	<b>84012</b>	50 ml	Sterile	Bulk	25	Blue Grad. Conical Test Tubes, 50 ml	<b>NEW</b>
	<b>84013</b>	50 ml	Sterile	Indiv. wrap	1	Blue Grad. Conical Test Tubes, 50 ml	<b>NEW</b>
	<b>84014</b>	50 ml	n/a	PS foam rack	25	Blue Grad. Conical Test Tubes, 50 ml	<b>NEW</b>
	<b>84015</b>	50 ml	Sterile	PS foam rack	25	Blue Grad. Conical Test Tubes, 50 ml	<b>NEW</b>