

# Kartell

## Labware division

---

Boxed Racks, Refill Racks and Stackable Racks  
for **Kartell** Tips



[www.kartell.com.au](http://www.kartell.com.au)



### RACKS - BOXED (Disposable)

Most KARTELL tips are available in autoclavable polypropylene racks which can be used with single pipettors but which are primarily designed for multi-channel units.

The configuration of these boxed racks which hold 96 tips, means that 4,8 and 12 channel multi-pipettors can easily be used.

	Art.	Tip Art.	Colour	Type	Qty of Tips	Packing	Notes
<b>CAPACITY 0.5 - 10 µL</b>	9646	956	neutral	Biohit®, Gilson® Micro	96	box & rack	Refill Avail.
<b>CAPACITY 0.5 - 10 µL</b>	9650	959	neutral	Nichiryo®, Eppendorf® Cristall, Kartell®, Socorex®, Gilson®	96	box & rack	Refill Avail.
<b>CAPACITY 2 - 20 µL</b>	9657	979	neutral	Gilson®, Eppendorf®	96	box & rack	Refill Avail.
<b>CAPACITY 2 - 200 µL</b>	9641	960	yellow	Universal	96	box & rack	Refill Avail.
<b>CAPACITY 2 - 200 µL</b>	9642	958	yellow	Gilson®	96	box & rack	Refill Avail.
<b>CAPACITY 5 - 200 µL</b>	9643	940	neutral grad.	Gilson®	96	box & rack	Refill Avail.
<b>CAPACITY 2 - 300 µL</b>	9644	971	neutral	Finnpipette®, Biohit®, Gilson®	96	box & rack	Refill Avail.
<b>CAPACITY 10 - 300 µL</b>	9647	935	neutral grad.	Biohit®, Eppendorf®, Gilson®	96	box & rack	Refill Avail.
<b>CAPACITY 100 - 1000 µL</b>	9640	961	blue	Universal	96	box & rack	Refill Avail.
<b>CAPACITY 100 - 1000 µL</b>	9645	941	blue	Gilson®	96	box & rack	Refill Avail.
<b>CAPACITY 1000 - 5000 µL</b>	9648	984	light blue	Kartell® Pluripet std and Pluripet II, Gilson®, Socorex®, Nichiryo®	96	box & rack	
<b>CAPACITY 1000 - 5000 µL</b>	9649	964	neutral	Kartell® Pluripet std, Eppendorf®, Biohit®, Oxford®	96	box & rack	
<b>CAPACITY 1000 - 5500 µL</b>	9847	942	neutral	Brand®, Finnpipette®, Treff®	96	box & rack	
<b>CAPACITY 1000 - 10000 µL</b>	9848	943	neutral	Socorex®, Gilson®, Finnpipette®	96	box & rack	

See  
Page  
2



## RACKS - REFILL

Majority of tips are also available in refillable racks. The 96 tips are held in position in a support unit on an 8x12 format which fits into the base unit. Once the tips have been used the support unit can be discarded and replaced with a new one complete with tips. Polycarbonate base units are also available.

	Art.	Tip Art.	Colour	Type	Qty of Tips	Packing
<b>CAPACITY 0.5 - 10 µL</b>	9747	956	neutral	Biohit®, Gilson® Micro	96	refill rack
<b>CAPACITY 0.5 - 10 µL</b>	9746	959	neutral	Nichiryo®, Eppendorf® Cristall, Kartell®, Socorex®, Gilson®	96	refill rack
<b>CAPACITY 2 - 20 µL</b>	9745	979	neutral	Gilson®, Eppendorf®	96	refill rack
<b>CAPACITY 2 - 200 µL</b>	9749	960	yellow	Universal	96	refill rack
<b>CAPACITY 2 - 200 µL</b>	9750	958	yellow	Gilson®	96	refill rack
<b>CAPACITY 5 - 200 µL</b>	9751	940	neutral grad.	Gilson®	96	refill rack
<b>CAPACITY 2 - 300 µL</b>	9752	971	neutral	Finnpipette®, Biohit®, Gilson®	96	refill rack
<b>CAPACITY 10 - 300 µL</b>	9754	935	neutral grad.	Biohit®, Eppendorf®, Gilson®	96	refill rack
<b>CAPACITY 100 - 1000 µL</b>	9748	961	blue	Universal	96	refill rack
<b>CAPACITY 100 - 1000 µL</b>	9753	941	blue	Gilson®	96	refill rack



## ECONOMY BOXED RACKS

Material: PC/PP

Economy option with one single Polycarbonate base and Polypropylene lid for each boxed rack.

	Art.	Tip Art.	Colour	Type	Qty of Tips	Packing
<b>CAPACITY 0.5 - 10 <math>\mu</math>L</b>	9624	959	crystal	Eppendorf® cristall	96	economy box & rack
<b>CAPACITY 2 - 200 <math>\mu</math>L</b>	9605	960	yellow	Universal	96	economy box & rack
<b>CAPACITY 2 - 200 <math>\mu</math>L</b>	9585	958	yellow	Gilson®	96	economy box & rack
<b>CAPACITY 200 - 1000 <math>\mu</math>L</b>	9615	957	blue	Universal	96	economy box & rack



## TIPS IN STACKED RACKS

Material: PC/PP

Space saving option with one single Polycarbonate base and Polypropylene lid for each stack. Made of 10 layers (2 - 200  $\mu$ L) or 5 layers (200 - 1000  $\mu$ L).

	Art.	Tip Art.	Colour	Type	Qty of Tips	Packing
<b>CAPACITY 2 - 200 <math>\mu</math>L</b>	9626	960	yellow	Universal	96 x 10 layers = 960 tips	stackable rack
<b>CAPACITY 2 - 200 <math>\mu</math>L</b>	9627	958	yellow	Gilson®	96 x 10 layers = 960 tips	stackable rack
<b>CAPACITY 200 - 1000 <math>\mu</math>L</b>	9628	957	blue	Universal	96 x 5 layers = 480 tips	stackable rack