



FIXED VOLUME MICRO PIPETTES - LABOFIX

In order to keep pace with the constantly increasing demand of fixed volume micropipettes with good quality standards, we are pleased to announce a new line called **LABOFIX**. Their structure completes the successful adjustable volume Pluripet II line and there are 10 models offering a wide choice of volume range from 5µl to 1.000µl.

Technical Features

- Completely autoclavable
- Suitable for the end-user in the field of Molecular Biology, Microbiology, Immunology, Cell culture, Analytical Chemistry, Biochemistry, Genetics etc. and for several applications
- CE Marked
- Accuracy and precision
- Light, ergonomic, built-in, streamlined tip ejector facilitates easy tip ejection and access to narrow necked bottles and tube.
- Constructed of fully autoclavable plastic and combined with stainless steel piston
- Individually Calibrated as per procedure laid down in ISO 8655 Standards;
- A calibration report is included with each pipette.

Art.	Volume	Accuracy	CV
		+ / - %	+ / - %
1369	5 µl	1.5	0.8
1370	10 µl	1.0	0.5
1371	20 µl	0.5	0.3
1372	25 µl	0.5	0.3
1373	50 µl	0.3	0.3
1374	100 µl	0.3	0.2
1375	200 µl	0.3	0.2
1376	250 µl	0.3	0.2
1377	500 µl	0.3	0.2
1378	1000 µl	0.3	0.2

	Art.	volume	Kartell tips code	
FIXED VOLUME MICRO PIPETTES - LABOFIX	1369	5 µl	956/959	NEW
	1370	10 µl	956/959	NEW
	1371	20 µl	958/960	NEW
	1372	25 µl	958/960	NEW
	1373	50 µl	958/960	NEW
	1374	100 µl	958/960	NEW
	1375	200 µl	958/960	NEW
	1376	250 µl	961	NEW
	1377	500 µl	961	NEW
	1378	1000 µl	961	NEW