



Organic chemistry catalogue

Chem-Supply has released the Acros Organics 2009 Price List and now stocks more products in Australia and offers faster indent lead-times.

Over 18,000 Acros compounds — available in 35,000 convenient pack sizes — can now be searched and ordered online at www.chemsupply.com.au.

Chem-Supply offers the full Acros product range of Essentials (for synthesis, separation and analysis), Building Blocks (used to construct more complex molecules) and Functional Reagents (to perform transformations of building blocks into target molecules).

Extra dry solvents and many other reagents come in the AcroSeal bottle, complete with triple sandwich septum — to ensure product integrity as well as safety and convenience. Also available is the range of Maybridge drug discovery compounds and screening libraries.

Chem-Supply Pty Ltd

Contact info and more items like this at wf.net.au/A516

Wash bottles

The Kartell Plastilab range of wash bottles includes traditional wash bottles, wide and narrow neck, graduated and ungraduated, capped and uncapped bottles as well as spare screw caps and adhesive labels which are sold separately.

The Kartell Traditional Wash Bottles, Article numbers: 179, 180, 182, 183 and 185, are moulded in a soft, yet resistant polyethylene to easily control liquid dispensing. The bottom of the bottle is designed to grant high stability during filling. The dispensing tip can be removed to increase flow. The sizes range from 50 to 1000 mL and are ungraduated. The Wide Mouth Wash Bottles, Article numbers: 1637 and 1638, are colour coded for easy identification and feature a wide mouth for safe filling. The dispensing tip is moulded in a single piece with the cap to prevent any leakage. The bottle's integral moulded tube allows easy control on fine tuning of liquid flux. These graduated wash bottles are available in 250 and 500 mL in 00 neutral, 04 blue, 06 yellow and 10 red. The Integral Wash Bottles, Article numbers: 1633 and 1634, are easy-squeeze wash bottles with an integral moulded tube which allows easy control on fine tuning of liquid flux. The dispensing tip is provided with a polyethylene closure cap. The ungraduated bottles are available in two sizes- 250 and 500 mL.

The range also includes Narrow Neck Bottles (without caps), that are suitable for use with water and oils - there are no shoulders to restrict the flow. The sizes range from 50 to 2000 mL and are graduated.

Spare screw caps and adhesive labels for the wash bottles are available.

Sieper & Co Pty Ltd

Contact info and more items like this at wf.net.au/A509





A09 A10

A11 A12



ASK NOW FOR THE LATEST BAC CATALOGUES FREE LAYOUT DESIGN WITH 3D MODELLING








BAC Storage Modules
WINNING 2008
INNOVATION AWARD
AUSTRALIA

BAC Modular Workplace and Storage Systems

(P) 02 9832 2777 • 193 - 195 Power Street Glendenning NSW 2761

www.bacsystems.com.au